

SIEMENS MAGNETOM Allegra syngo MR 2004A

\\USER\Narr\Depression\Brain\Circle Scout

Scan Time: 0:17 Voxel size: 1.2x1.2x7.0 [mm] Rel. SNR: 1.00 SIEMENS: gre_circle

Routine

Slice group 1	
Slices	1
Dist. factor	20 [%]
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0 [deg]
Slice group 2	
Slices	1
Dist. factor	20 [%]
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Rotation	0 [deg]
Slice group 3	
Slices	1
Dist. factor	20 [%]
Position	Isocenter
Orientation	Coronal
Phase enc. dir.	R >> L
Rotation	0 [deg]
Phase oversampling	0 [%]
FoV read	300 [mm]
FoV phase	100.0 [%]
Slice thickness	7 [mm]
TR	20 [ms]
TE	3.39 [ms]
Averages	1
Concatenations	3
Filter	None
Coil elements	HE

Contrast

TD	0 [ms]
MTC	0
Magn. preparation	None
Flip angle	25 [deg]
Reconstruction	Magnitude
Fat suppr.	None
Water suppr.	None
Measurements	1

Resolution

Base resolution	256
Phase resolution	100 [%]
Phase partial Fourier	Off
Filter 1	
Raw filter	Off
Filter 2	
Large FoV	Off
Filter 3	
Normalize	Off
Filter 4	
Elliptical filter	Off
Interpolation	0

PAT mode	None

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Saturation mode	Standard
Special sat.	None

System

Save uncombined	0
Scan at current TP	1
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Head 3T / HE	1

Shim mode	Tune up
Confirm freq. adjustment	0
Assume Silicone	0
Ref. amplitude [1H]	140.000 [V]
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0 [deg]
R >> L	350 [mm]
A >> P	263 [mm]
F >> H	350 [mm]

Physio

1st Signal/Mode	None
Segments	1

Dark blood	0

Inline

Subtract	0
Std-Dev-Sag	0
Std-Dev-Cor	0
Std-Dev-Tra	0
Std-Dev-Time	0
MIP-Sag	0
MIP-Cor	0
MIP-Tra	0
MIP-Time	0
Save original images	1

Sequence

Introduction	1
Dimension	2D
Phase stabilisation	0
Averaging mode	Short term
Asymmetric echo	Off
Contrasts	1
Bandwidth	260 [Hz/Px]
Flow comp.	No

RF pulse type	Normal
Gradient mode	Fast
Excitation	Slice-sel.
RF spoiling	1

SIEMENS MAGNETOM Allegra syngo MR 2004A

\\USER\Narr\Depression\Brain\tfl_mgh_me_vNav

Scan Time: 8:07 Voxel size: 1.3x1.0x1.3 [mm] Rel. SNR: 1.00 USER: epi_nav\tfl_mgh_multiecho_epinav_r5

Routine

Slab group 1	
Slabs	1
Dist. factor	50 [%]
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Rotation	0 [deg]
Phase oversampling	0 [%]
Slice oversampling	0 [%]
Slices per slab	128
FoV read	256 [mm]
FoV phase	100.0 [%]
Slice thickness	1.33 [mm]
TR	2530 [ms]
TE[1]	1.53 [ms]
TE[2]	3.21 [ms]
TE[3]	4.89 [ms]
TE[4]	6.57 [ms]
Averages	1
Concatenations	1
Filter	None
Coil elements	HE

Contrast

Magn. preparation	Non-sel. IR
T1	1100 [ms]
Flip angle	7 [deg]
Reconstruction	Magnitude
Fat suppr.	None
Water suppr.	None
Measurements	1

Resolution

Base resolution	256
Phase resolution	75 [%]
Slice resolution	100 [%]
Phase partial Fourier	Off
Slice partial Fourier	Off
Filter 1	
Raw filter	Off
Filter 2	
Large FoV	Off
Filter 3	
Normalize	Off
Filter 4	
Elliptical filter	Off
Interpolation	0

PAT mode	None

Geometry

Multi-slice mode	Single shot
Series	Interleaved

System

Save uncombined	0
Scan at current TP	0
Scan region position	H
Scan region position	0 [mm]
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Head 3T / HE	1

Shim mode	Standard

Confirm freq. adjustment	0
Assume Silicone	0
Ref. amplitude [1H]	140.000 [V]
Adjust volume	
Position	Isocenter
Orientation	Sagittal
Rotation	0 [deg]
F >> H	256 [mm]
A >> P	256 [mm]
R >> L	171 [mm]

Physio

1st Signal/Mode	None

Dark blood	0

Inline

Subtract	0
Std-Dev-Sag	0
Std-Dev-Cor	0
Std-Dev-Tra	0
Std-Dev-Time	0
MIP-Sag	0
MIP-Cor	0
MIP-Tra	0
MIP-Time	0
Save original images	1

Sequence

Introduction	1
Dimension	3D
Elliptical scanning	0
Averaging mode	Long term
Asymmetric echo	Off
Contrasts	4
Bandwidth	651 [Hz/Px]
Flow comp.	No
Echo spacing	8.9 [ms]

RF pulse type	Fast
Gradient mode	Fast
Excitation	Non-sel.
RF spoiling	1

Add. scale factor	4
Gradient spoiling	Siemens
Gradient moment factor	1
Siemens reconstruction	1
Save raw k-space data	0
Averaging	RMS

SIEMENS MAGNETOM Allegra syngo MR 2004A

\\USER\Narr\Depression\Brain\Matched Bandwidth Hi-Res

Scan Time: 1:30 Voxel size: 1.7x1.7x4.0 [mm] Rel. SNR: 1.00 SIEMENS: ep_seg_se

Routine

Slice group 1	
Slices	34
Dist. factor	25 [%]
Position	R1.4 P6.8 H10.4 [mm]
Orientation	T > C-7.0
Phase enc. dir.	A >> P
Rotation	0 [deg]
Phase oversampling	0 [%]
FoV read	220 [mm]
FoV phase	100.0 [%]
Slice thickness	4 [mm]
TR	5000 [ms]
TE	33 [ms]
Averages	4
Concatenations	1
Filter	None
Coil elements	HE

R >> L	220 [mm]
A >> P	220 [mm]
F >> H	169 [mm]

Physio

1st Signal/Mode	None
Resp. control	Off

Sequence

Introduction	0
Dimension	2D
Averaging mode	Long term
Bandwidth	1302 [Hz/Px]
Free echo spacing	0
Echo spacing	0.85 [ms]
EPI factor	33
RF pulse type	Normal
Gradient mode	Fast

Contrast

MTC	0
Magn. preparation	None
Flip angle	90 [deg]
Reconstruction	Magnitude
Fat suppr.	Fat sat.
Measurements	1

Resolution

Base resolution	128
Phase resolution	100 [%]
Phase partial Fourier	Off
Filter 1	
Raw filter	Off
Filter 2	
Large FoV	Off
Filter 3	
Normalize	Off
Filter 4	
Elliptical filter	Off
Interpolation	0

Geometry

Multi-slice mode	Interleaved
Series	Interleaved
Special sat.	None

System

Save uncombined	0
Scan at current TP	0
Scan region position	H
Scan region position	0 [mm]
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Head 3T / HE	1
Shim mode	Standard
Confirm freq. adjustment	0
Assume Silicone	0
Ref. amplitude [1H]	140.000 [V]
Adjust volume	
Position	R1.4 P6.8 H10.4 [mm]
Orientation	T > C-7.0
Rotation	0 [deg]

SIEMENS MAGNETOM Allegra syngo MR 2004A

\\USER\Narr\Depression\Brain\Resting

Scan Time: 6:04 Voxel size: 3.4x3.4x4.0 [mm] Rel. SNR: 1.00 SIEMENS: ep2d_bold

Routine	1st Signal/Mode	None
Slice group 1		
Slices	34	
Dist. factor	25 [%]	
Position	R1.4 P6.8 H10.4 [mm]	
Orientation	T > C-7.0	
Phase enc. dir.	A >> P	
Rotation	0 [deg]	
Phase oversampling	0 [%]	
FoV read	220 [mm]	
FoV phase	100.0 [%]	
Slice thickness	4 [mm]	
TR	2000 [ms]	
TE	30 [ms]	
Averages	1	
Concatenations	1	
Filter	None	
Coil elements	HE	
Contrast		
MTC	0	
Flip angle	77 [deg]	
Reconstruction	Magnitude	
Fat suppr.	Fat sat.	
Measurements	180	
Delay in TR	0 [ms]	
Multiple series	0	
Resolution		
Base resolution	64	
Phase resolution	100 [%]	
Phase partial Fourier	Off	
Filter 1		
Raw filter	Off	
Interpolation	0	

PAT mode	None	
Geometry		
Multi-slice mode	Interleaved	
Series	Interleaved	

Special sat.	None	
System		
Scan at current TP	0	
Scan region position	H	
Scan region position	0 [mm]	
MSMA	S - C - T	
Sagittal	R >> L	
Coronal	A >> P	
Transversal	F >> H	
Head 3T / HE	1	

Shim mode	Standard	
Confirm freq. adjustment	0	
Assume Silicone	0	
Ref. amplitude [1H]	140.000 [V]	
Adjust volume		
Position	R1.4 P6.8 H10.4 [mm]	
Orientation	T > C-7.0	
Rotation	0 [deg]	
R >> L	220 [mm]	
A >> P	220 [mm]	
F >> H	169 [mm]	
Physio		

BOLD		
t-Test	0	
Threshold	3.00	
Window	Growing	
Dynamic t-maps	1	
Starting ignore meas	0	
Paradigm size	12	
Meas[1]	Baseline	
Meas[2]	Active	
Meas[3]	Baseline	
Meas[4]	Active	
Meas[5]	Baseline	
Meas[6]	Active	
Meas[7]	Baseline	
Meas[8]	Active	
Meas[9]	Baseline	
Meas[10]	Active	
Meas[11]	Baseline	
Meas[12]	Active	
Motion correction	0	
Spatial filter	0	
Sequence		
Introduction	0	
Averaging mode	Long term	
Bandwidth	2604 [Hz/Px]	
Free echo spacing	0	
Echo spacing	0.44 [ms]	

EPI factor	64	
RF pulse type	Normal	
Gradient mode	Fast	

SIEMENS MAGNETOM Allegra syngo MR 2004A

\\USER\Narr\Depression\Brain\ep2d_casl_UI

Scan Time: 4:08 Voxel size: 4.0x4.0x6.0 [mm] Rel. SNR: 1.00 USER: ep2d_casl_UI_iPAT2

Routine

Slice group 1	
Slices	18
Dist. factor	25 [%]
Position	L2.9 A3.1 H29.3 [mm]
Orientation	T > C13.0
Phase enc. dir.	A >> P
Rotation	0 [deg]
Phase oversampling	0 [%]
FoV read	256 [mm]
FoV phase	100.0 [%]
Slice thickness	6 [mm]
TR	4000 [ms]
TE	16 [ms]
Averages	1
Concatenations	1
Filter	None
Coil elements	HE

Contrast

MTC	0
Flip angle	90 [deg]
Reconstruction	Magnitude
Fat suppr.	Fat sat.
Measurements	60
Delay in TR	0 [ms]
Multiple series	0

Resolution

Base resolution	64
Phase resolution	100 [%]
Phase partial Fourier	Off
Filter 1	
Raw filter	Off
Interpolation	0

PAT mode	None

Geometry

Multi-slice mode	Interleaved
Series	Ascending

Special sat.	None

System

Scan at current TP	0
Scan region position	H
Scan region position	0 [mm]
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Head 3T / HE	1

Shim mode	Standard
Confirm freq. adjustment	0
Assume Silicone	0
Ref. amplitude [1H]	140.000 [V]
Adjust volume	
Position	L2.9 A3.1 H29.3 [mm]
Orientation	T > C13.0
Rotation	0 [deg]
R >> L	256 [mm]
A >> P	256 [mm]
F >> H	134 [mm]

Physio

1st Signal/Mode

None

BOLD

t-Test	0
Threshold	4.00
Window	Growing
Dynamic t-maps	0
Starting ignore meas	0
Paradigm size	20
Meas[1]	Ignore
Meas[2]	Ignore
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Ignore
Meas[12]	Ignore
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	1
Interpolation	3D-K-space
Spatial filter	0

Sequence

Introduction	0
Averaging mode	Long term
Bandwidth	3004 [Hz/Px]
Free echo spacing	0
Echo spacing	0.38 [ms]

EPI factor	64
RF pulse type	Normal
Gradient mode	Fast

CASL Method	Multi-slice
Label Offset	80 mm
Post Label Delay	1000000 us
Num RF Blocks	20
RF GAP	5000 us
Pre Label Delay	0 us

SIEMENS MAGNETOM Allegra syngo MR 2004A

\\USER\Narr\Depression\Brain\ep_set_shmocoSVS_5mm_adult

+ Scan Time: 1.2 [s] Voxel size: 5.0x5.0x5.0 [mm] Rel. SNR: 1.00 USER: epi_nav\lep_moco_set

Routine

Slab group 1	
Slabs	1
Dist. factor	50 [%]
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	R >> L
Rotation	90 [deg]
Phase oversampling	0 [%]
Slice oversampling	0 [%]
Slices per slab	28
FoV read	220 [mm]
FoV phase	90.9 [%]
Slice thickness	5 [mm]
TR	21 [ms]
TE[1]	8 [ms]
TE[2]	12.8 [ms]
Averages	1
Concatenations	1
Filter	None
Coil elements	HE

Rotation	0 [deg]
R >> L	350 [mm]
A >> P	263 [mm]
F >> H	350 [mm]

Physio

1st Signal/Mode	None
-----------------	------

Sequence

Introduction	0
Dimension	3D
Averaging mode	Long term
Contrasts	2
Bandwidth	3918 [Hz/Px]
Free echo spacing	0
Echo spacing	0.31 [ms]

EPI factor	40
RF pulse type	Normal
Gradient mode	Normal
RF spoiling	1

Shim Navigator	1
Delta TE	4800
Protocol Name	SVS

Contrast

MTC	0
Flip angle	2 [deg]
Reconstruction	Magn./Phase
Fat suppr.	None
Measurements	1

Resolution

Base resolution	44
Phase resolution	100 [%]
Slice resolution	100 [%]
Phase partial Fourier	Off
Slice partial Fourier	Off
Filter 1	
Raw filter	Off
Filter 2	
Large FoV	Off
Filter 3	
Normalize	Off
Filter 4	
Elliptical filter	Off
Interpolation	0

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat.	None

System

Save uncombined	0
Scan at current TP	1
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Head 3T / HE	1

Shim mode	Tune up
Confirm freq. adjustment	0
Assume Silicone	0
Ref. amplitude [1H]	140.000 [V]
Adjust volume	
Position	Isocenter
Orientation	Transversal

SIEMENS MAGNETOM Allegra syngo MR 2004A

\\USER\Narr\Depression\Brain\svs_vNav_WS Lt Hippo

Scan Time: 4:50 Vol: 30 x12 x12 [mm] Rel. SNR: 1.00 USER: epi_nav\svs_se_NAV_VA25_r5

Routine

Position	L28.0 P2.5 F41.2 [mm]
Orientation	T > C21.0
Rotation	0 [deg]
Vol A >> P	30 [mm]
Vol R >> L	12 [mm]
Vol F >> H	12 [mm]
TR	2200 [ms]
TE	30 [ms]
Averages	128
Coil elements	HE

Contrast

Flip angle	90 [deg]
Water suppr.	Water sat.
Water suppr. BW	70 [Hz]
Measurements	1

Resolution

Vector size	1024
-------------	------

Geometry

System

Save uncombined	1
Scan at current TP	0
Scan region position	H
Scan region position	0 [mm]
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Head 3T / HE	1

Shim mode	Advanced
Adj. water suppr.	0
Confirm freq. adjustment	0
Assume Silicone	0
Ref. amplitude [1H]	140.000 [V]
Adjust volume	
Position	L28.0 P2.5 F41.2 [mm]
Orientation	T > C21.0
Rotation	90 [deg]
A >> P	30 [mm]
R >> L	12 [mm]
F >> H	12 [mm]

Sequence

Preparation scans	4
Delta frequency	-3.4 [ppm]
Phase cycling	Auto
Bandwidth	2000 [Hz]
Acquisition duration	512 [ms]

EPI Navigator On	1
Apply Feedback	SVS Shim
Feedback Wait	170000
Wait Remain	10000
.	0
Empty Ice	0
Disable Nav Shim	0
Nav Prot	user_SVS

SIEMENS MAGNETOM Allegra syngo MR 2004A

\\USER\Narr\Depression\Brain\svs_vNav_NWS Lt Hippo

Scan Time: 0:32 Vol: 30 x12 x12 [mm] Rel. SNR: 1.00 USER: epi_nav\svs_se_NAV_VA25_r5

Routine

Position	L28.0 P2.5 F41.2 [mm]
Orientation	T > C21.0
Rotation	0 [deg]
Vol A >> P	30 [mm]
Vol R >> L	12 [mm]
Vol F >> H	12 [mm]
TR	8000 [ms]
TE	30 [ms]
Averages	1
Coil elements	HE

Contrast

Flip angle	90 [deg]
Water suppr.	None
Water suppr. BW	35 [Hz]
Measurements	1

Resolution

Vector size	1024
-------------	------

Geometry

System

Save uncombined	1
Scan at current TP	0
Scan region position	H
Scan region position	0 [mm]
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Head 3T / HE	1

Shim mode	Advanced
Adj. water suppr.	0
Confirm freq. adjustment	0
Assume Silicone	0
Ref. amplitude [1H]	140.000 [V]
Adjust volume	
Position	L28.0 P2.5 F41.2 [mm]
Orientation	T > C21.0
Rotation	90 [deg]
A >> P	30 [mm]
R >> L	12 [mm]
F >> H	12 [mm]

Sequence

Preparation scans	3
Delta frequency	-3.4 [ppm]
Phase cycling	Auto
Bandwidth	2000 [Hz]
Acquisition duration	512 [ms]

EPI Navigator On	1
Apply Feedback	SVS Shim
Feedback Wait	170000
Wait Remain	10000
.	0
Empty Ice	0
Disable Nav Shim	0
Nav Prot	user_SVS

SIEMENS MAGNETOM Allegra syngo MR 2004A

\\USER\Nar\Depression\Brain\svs_vNav_WS Rt Hippo

Scan Time: 4:50 Vol: 30 x12 x12 [mm] Rel. SNR: 1.00 USER: epi_nav\svs_se_NAV_VA25_r5

Routine

Position	R20.0 P26.0 F39.0 [mm]
Orientation	T > C16.6
Rotation	0 [deg]
Vol A >> P	30 [mm]
Vol R >> L	12 [mm]
Vol F >> H	12 [mm]
TR	2200 [ms]
TE	30 [ms]
Averages	128
Coil elements	HE

Contrast

Flip angle	90 [deg]
Water suppr.	Water sat.
Water suppr. BW	70 [Hz]
Measurements	1

Resolution

Vector size	1024
-------------	------

Geometry

System

Save uncombined	1
Scan at current TP	0
Scan region position	H
Scan region position	0 [mm]
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Head 3T / HE	1

Shim mode	Advanced
Adj. water suppr.	0
Confirm freq. adjustment	0
Assume Silicone	0
Ref. amplitude [1H]	140.000 [V]
Adjust volume	
Position	R20.0 P26.0 F39.0 [mm]
Orientation	T > C16.6
Rotation	90 [deg]
A >> P	30 [mm]
R >> L	12 [mm]
F >> H	12 [mm]

Sequence

Preparation scans	4
Delta frequency	-3.4 [ppm]
Phase cycling	Auto
Bandwidth	2000 [Hz]
Acquisition duration	512 [ms]

EPI Navigator On	1
Apply Feedback	SVS Shim
Feedback Wait	170000
Wait Remain	10000
.	0
Empty Ice	0
Disable Nav Shim	0
Nav Prot	user_SVS

SIEMENS MAGNETOM Allegra syngo MR 2004A

\\USER\Narr\Depression\Brain\svs_vNav_NWS Rt Hippo

Scan Time: 0:32 Vol: 30 x12 x12 [mm] Rel. SNR: 1.00 USER: epi_nav\svs_se_NAV_VA25_r5

Routine

Position	R20.0 P26.0 F39.0 [mm]
Orientation	T > C16.6
Rotation	0 [deg]
Vol A >> P	30 [mm]
Vol R >> L	12 [mm]
Vol F >> H	12 [mm]
TR	8000 [ms]
TE	30 [ms]
Averages	1
Coil elements	HE

Contrast

Flip angle	90 [deg]
Water suppr.	None
Water suppr. BW	35 [Hz]
Measurements	1

Resolution

Vector size	1024
-------------	------

Geometry

System

Save uncombined	1
Scan at current TP	0
Scan region position	H
Scan region position	0 [mm]
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Head 3T / HE	1

Shim mode	Advanced
Adj. water suppr.	0
Confirm freq. adjustment	0
Assume Silicone	0
Ref. amplitude [1H]	140.000 [V]
Adjust volume	
Position	R20.0 P26.0 F39.0 [mm]
Orientation	T > C16.6
Rotation	90 [deg]
A >> P	30 [mm]
R >> L	12 [mm]
F >> H	12 [mm]

Sequence

Preparation scans	3
Delta frequency	-3.4 [ppm]
Phase cycling	Auto
Bandwidth	2000 [Hz]
Acquisition duration	512 [ms]

EPI Navigator On	1
Apply Feedback	SVS Shim
Feedback Wait	170000
Wait Remain	10000
.	0
Empty Ice	0
Disable Nav Shim	0
Nav Prot	user_SVS

SIEMENS MAGNETOM Allegra syngo MR 2004A

\\USER\Narr\Depression\Brain\svs_vNav_WS DorCing

Scan Time: 4:50 Vol: 20 x18 x12 [mm] Rel. SNR: 1.00 USER: epi_nav\svs_se_NAV_VA25_r5

Routine

Position	R14.0 A45.0 F3.5 [mm]
Orientation	T > C29.9
Rotation	0 [deg]
Vol A >> P	20 [mm]
Vol R >> L	18 [mm]
Vol F >> H	12 [mm]
TR	2200 [ms]
TE	30 [ms]
Averages	128
Coil elements	HE

Contrast

Flip angle	90 [deg]
Water suppr.	Water sat.
Water suppr. BW	70 [Hz]
Measurements	1

Resolution

Vector size	1024
-------------	------

Geometry

System

Save uncombined	1
Scan at current TP	0
Scan region position	H
Scan region position	0 [mm]
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Head 3T / HE	1

Shim mode	Advanced
Adj. water suppr.	0
Confirm freq. adjustment	0
Assume Silicone	0
Ref. amplitude [1H]	140.000 [V]
Adjust volume	
Position	R14.0 A45.0 F3.5 [mm]
Orientation	T > C29.9
Rotation	90 [deg]
A >> P	20 [mm]
R >> L	18 [mm]
F >> H	12 [mm]

Sequence

Preparation scans	4
Delta frequency	-3.4 [ppm]
Phase cycling	Auto
Bandwidth	2000 [Hz]
Acquisition duration	512 [ms]

EPI Navigator On	1
Apply Feedback	SVS Shim
Feedback Wait	170000
Wait Remain	10000
.	0
Empty Ice	0
Disable Nav Shim	0
Nav Prot	user_SVS

SIEMENS MAGNETOM Allegra syngo MR 2004A

\\USER\Narr\Depression\Brain\svs_vNav_NWS DorCing

Scan Time: 0:32 Vol: 20 x18 x12 [mm] Rel. SNR: 1.00 USER: epi_nav\svs_se_NAV_VA25_r5

Routine

Position	R14.0 A45.0 F3.5 [mm]
Orientation	T > C29.9
Rotation	0 [deg]
Vol A >> P	20 [mm]
Vol R >> L	18 [mm]
Vol F >> H	12 [mm]
TR	8000 [ms]
TE	30 [ms]
Averages	1
Coil elements	HE

Contrast

Flip angle	90 [deg]
Water suppr.	None
Water suppr. BW	35 [Hz]
Measurements	1

Resolution

Vector size	1024
-------------	------

Geometry

System

Save uncombined	1
Scan at current TP	0
Scan region position	H
Scan region position	0 [mm]
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Head 3T / HE	1

Shim mode	Advanced
Adj. water suppr.	0
Confirm freq. adjustment	0
Assume Silicone	0
Ref. amplitude [1H]	140.000 [V]
Adjust volume	
Position	R14.0 A45.0 F3.5 [mm]
Orientation	T > C29.9
Rotation	90 [deg]
A >> P	20 [mm]
R >> L	18 [mm]
F >> H	12 [mm]

Sequence

Preparation scans	3
Delta frequency	-3.4 [ppm]
Phase cycling	Auto
Bandwidth	2000 [Hz]
Acquisition duration	512 [ms]

EPI Navigator On	1
Apply Feedback	SVS Shim
Feedback Wait	170000
Wait Remain	10000
.	0
Empty Ice	0
Disable Nav Shim	0
Nav Prot	user_SVS

SIEMENS MAGNETOM Allegra syngo MR 2004A

\\USER\Narr\Depression\Brain\svs_vNav_WS subCing

Scan Time: 4:50 Vol: 20 x18 x12 [mm] Rel. SNR: 1.00 USER: epi_nav\svs_se_NAV_VA25_r5

Routine

Position	R4.0 A51.0 F35.0 [mm]
Orientation	T > C-41.0
Rotation	0 [deg]
Vol A >> P	20 [mm]
Vol R >> L	18 [mm]
Vol F >> H	12 [mm]
TR	2200 [ms]
TE	30 [ms]
Averages	128
Coil elements	HE

Contrast

Flip angle	90 [deg]
Water suppr.	Water sat.
Water suppr. BW	70 [Hz]
Measurements	1

Resolution

Vector size	1024
-------------	------

Geometry

System

Save uncombined	1
Scan at current TP	0
Scan region position	H
Scan region position	0 [mm]
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Head 3T / HE	1

Shim mode	Advanced
Adj. water suppr.	0
Confirm freq. adjustment	0
Assume Silicone	0
Ref. amplitude [1H]	140.000 [V]
Adjust volume	
Position	R4.0 A51.0 F35.0 [mm]
Orientation	T > C-41.0
Rotation	90 [deg]
A >> P	20 [mm]
R >> L	18 [mm]
F >> H	12 [mm]

Sequence

Preparation scans	4
Delta frequency	-3.4 [ppm]
Phase cycling	Auto
Bandwidth	2000 [Hz]
Acquisition duration	512 [ms]

EPI Navigator On	1
Apply Feedback	SVS Shim
Feedback Wait	170000
Wait Remain	10000
.	0
Empty Ice	0
Disable Nav Shim	0
Nav Prot	user_SVS

SIEMENS MAGNETOM Allegra syngo MR 2004A

\\USER\Narr\Depression\Brain\svs_vNav_NWS subCing

Scan Time: 0:32 Vol: 20 x18 x12 [mm] Rel. SNR: 1.00 USER: epi_nav\svs_se_NAV_VA25_r5

Routine

Position	R4.0 A51.0 F35.0 [mm]
Orientation	T > C-41.0
Rotation	0 [deg]
Vol A >> P	20 [mm]
Vol R >> L	18 [mm]
Vol F >> H	12 [mm]
TR	8000 [ms]
TE	30 [ms]
Averages	1
Coil elements	HE

Contrast

Flip angle	90 [deg]
Water suppr.	None
Water suppr. BW	35 [Hz]
Measurements	1

Resolution

Vector size	1024
-------------	------

Geometry

System

Save uncombined	1
Scan at current TP	0
Scan region position	H
Scan region position	0 [mm]
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Head 3T / HE	1

Shim mode	Advanced
Adj. water suppr.	0
Confirm freq. adjustment	0
Assume Silicone	0
Ref. amplitude [1H]	140.000 [V]
Adjust volume	
Position	R4.0 A51.0 F35.0 [mm]
Orientation	T > C-41.0
Rotation	90 [deg]
A >> P	20 [mm]
R >> L	18 [mm]
F >> H	12 [mm]

Sequence

Preparation scans	3
Delta frequency	-3.4 [ppm]
Phase cycling	Auto
Bandwidth	2000 [Hz]
Acquisition duration	512 [ms]

EPI Navigator On	1
Apply Feedback	SVS Shim
Feedback Wait	170000
Wait Remain	10000
.	0
Empty Ice	0
Disable Nav Shim	0
Nav Prot	user_SVS

SIEMENS MAGNETOM Allegra syngo MR 2004A

\\USER\Narr\Depression\Brain\ep2d_diff_WIP_96x96_37

+ Scan Time: 4:17 Voxel size: 2.5x2.5x2.5 [mm] Rel. SNR: 1.00 USER: ep2d_diff

Routine

Slice group 1	
Slices	55
Dist. factor	0 [%]
Position	L0.9 P5.7 H1.6 [mm]
Orientation	T > C-3.1 > S0.2
Phase enc. dir.	A >> P
Rotation	0 [deg]
Phase oversampling	0 [%]
FoV read	240 [mm]
FoV phase	100.0 [%]
Slice thickness	2.5 [mm]
TR	7300 [ms]
TE	95 [ms]
Averages	1
Concatenations	1
Filter	None
Coil elements	HE

Contrast

MTC	0
Magn. preparation	None
Reconstruction	Magnitude
Fat suppr.	Fat sat.
Measurements	1
Delay in TR	0 [ms]

Resolution

Base resolution	96
Phase resolution	100 [%]
Phase partial Fourier	Off
Filter 1	
Raw filter	Off
Filter 2	
Large FoV	Off
Filter 3	
Normalize	Off
Filter 4	
Elliptical filter	Off
Interpolation	1

PAT mode	None

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat.	None

System

Scan at current TP	0
Scan region position	H
Scan region position	0 [mm]
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Head 3T / HE	1

Shim mode	Standard
Confirm freq. adjustment	0
Assume Silicone	0
Ref. amplitude [1H]	140.000 [V]
Adjust volume	
Position	L0.9 P5.7 H1.6 [mm]
Orientation	T > C-3.1 > S0.2
Rotation	0 [deg]

R >> L	240 [mm]
A >> P	240 [mm]
F >> H	138 [mm]

Physio

1st Signal/Mode	None
-----------------	------

Diff

Diffusion mode	Free
Diff. weightings	1
b-value	1000 [s/mm ²]
Diff. weighted images	1
Trace weighted images	0
Average ADC maps	0
Individual ADC maps	0
Noise level	40
Diff. directions	34

Sequence

Introduction	1
Averaging mode	Long term
Bandwidth	2894 [Hz/Px]
Free echo spacing	0
Echo spacing	0.47 [ms]

EPI factor	96
RF pulse type	Normal
Gradient mode	Fast*

SIEMENS MAGNETOM Allegra syngo MR 2004A

\\USER\Narr\Depression\Brain\ep2d_diff_WIP_96x96_34

+ Scan Time: 4:17 Voxel size: 2.5x2.5x2.5 [mm] Rel. SNR: 1.00 USER: ep2d_diff

Routine

Slice group 1	
Slices	55
Dist. factor	0 [%]
Position	L0.9 P5.7 H1.6 [mm]
Orientation	T > C-3.1 > S0.2
Phase enc. dir.	A >> P
Rotation	0 [deg]
Phase oversampling	0 [%]
FoV read	240 [mm]
FoV phase	100.0 [%]
Slice thickness	2.5 [mm]
TR	7300 [ms]
TE	95 [ms]
Averages	1
Concatenations	1
Filter	None
Coil elements	HE

R >> L	240 [mm]
A >> P	240 [mm]
F >> H	138 [mm]

Physio

1st Signal/Mode	None
-----------------	------

Diff

Diffusion mode	Free
Diff. weightings	1
b-value	1000 [s/mm ²]
Diff. weighted images	1
Trace weighted images	0
Average ADC maps	0
Individual ADC maps	0
Noise level	40
Diff. directions	34

Sequence

Introduction	1
Averaging mode	Long term
Bandwidth	2894 [Hz/Px]
Free echo spacing	0
Echo spacing	0.47 [ms]

EPI factor	96
RF pulse type	Normal
Gradient mode	Fast*

Contrast

MTC	0
Magn. preparation	None
Reconstruction	Magnitude
Fat suppr.	Fat sat.
Measurements	1
Delay in TR	0 [ms]

Resolution

Base resolution	96
Phase resolution	100 [%]
Phase partial Fourier	Off
Filter 1	
Raw filter	Off
Filter 2	
Large FoV	Off
Filter 3	
Normalize	Off
Filter 4	
Elliptical filter	Off
Interpolation	1

PAT mode	None

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat.	None

System

Scan at current TP	0
Scan region position	H
Scan region position	0 [mm]
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Head 3T / HE	1

Shim mode	Standard
Confirm freq. adjustment	0
Assume Silicone	0
Ref. amplitude [1H]	140.000 [V]
Adjust volume	
Position	L0.9 P5.7 H1.6 [mm]
Orientation	T > C-3.1 > S0.2
Rotation	0 [deg]

SIEMENS MAGNETOM Allegra syngo MR 2004A

\\USER\Nar\Depression\Brain\gre_field_map

Scan Time: 4:20 Voxel size: 2.0x2.0x2.0 [mm] Rel. SNR: 1.00 SIEMENS: gre_field_mapping

Routine

Slice group 1	
Slices	54
Dist. factor	50 [%]
Position	R1.4 P6.1 H6.8 [mm]
Orientation	T > C-7.0
Phase enc. dir.	A >> P
Rotation	0 [deg]
Phase oversampling	0 [%]
FoV read	256 [mm]
FoV phase	100.0 [%]
Slice thickness	2 [mm]
TR	1000 [ms]
TE[1]	5.53 [ms]
TE[2]	7.99 [ms]
Averages	1
Concatenations	1
Filter	None
Coil elements	HE

A >> P 256 [mm]
F >> H 161 [mm]

Sequence

Introduction	1
Dimension	2D
Averaging mode	Short term
Asymmetric echo	Off
Contrasts	2
Bandwidth	1302 [Hz/Px]
Flow comp.	Yes

RF pulse type	Normal
Gradient mode	Fast
RF spoiling	1

Contrast

MTC	0
Flip angle	55 [deg]
Reconstruction	Magn./Phase
Fat suppr.	None
Measurements	1

Resolution

Base resolution	128
Phase resolution	100 [%]
Phase partial Fourier	Off
Filter 1	
Raw filter	Off
Filter 2	
Large FoV	Off
Filter 3	
Normalize	Off
Filter 4	
Elliptical filter	Off
Interpolation	0

Geometry

Multi-slice mode	Interleaved
Series	Interleaved
<hr style="border-top: 1px dashed black;"/>	
Special sat.	None

System

Save uncombined	0
Scan at current TP	0
Scan region position	H
Scan region position	0 [mm]
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Head 3T / HE	1
<hr style="border-top: 1px dashed black;"/>	
Shim mode	Standard
Confirm freq. adjustment	0
Assume Silicone	0
Ref. amplitude [1H]	140.000 [V]
Adjust volume	
Position	R1.4 P6.1 H6.8 [mm]
Orientation	T > C-7.0
Rotation	0 [deg]
R >> L	256 [mm]

Table of contents

\\USER	Narr	Depression	Brain
			Circle Scout
			tfl_mgh_me_vNav
			Matched Bandwidth Hi-Res
			Resting
			ep2d_casl_UI
			ep_set_shmocoSVS_5mm_adult
			svs_vNav_WS Lt Hippo
			svs_vNav_NWS Lt Hippo
			svs_vNav_WS Rt Hippo
			svs_vNav_NWS Rt Hippo
			svs_vNav_WS DorCing
			svs_vNav_NWS DorCing
			svs_vNav_WS subCing
			svs_vNav_NWS subCing
			ep2d_diff_WIP_96x96_37
			ep2d_diff_WIP_96x96_34
			gre_field_map